

Rehydrating Lyophilized Materials

Purpose The purpose of this protocol is to rehydrate the lyophilized (dried) biologicals and prepare stock solutions of **bacteria** and **control phage** to be stored and used at your institution. It is critical that cultures are prepared from a single colony and lysates prepared from a single plaque.

Supplies

- Lyophilized Vial
- Micropipette
- Pipette tips
- O-ring tube (-80°C compatible)
- Liquid Growth Media
- Agar Plates
- Sterile Loop
- 40% Glycerol
- DMSO
- Freezer ready labels (such as babies)

Getting Started

Prepare your bench for aseptic work and assemble your supplies.




Figure A: Photos of the lyophilized material. Preparing the vial **1**. The vial as it is received. **2**. The tab after it is lifted. **3**. Removal of the pull-tab from the vial. **4**. The freeze-dried, vacuum-sealed vial. **5**. The vial, prepared to receive media.

Procedures: Bacteria

A. Lyophilized Bacteria Rehydration

1. Remove packaging from vial (Fig. A-1).
2. Remove lid from vial (Fig. A-1 to -5).
3. With aseptic technique, transfer 500 μ L of liquid growth media into the vial (Fig. A-5).
4. The lyophilized cake will disintegrate almost immediately.
5. Gently pipette up and down 3-5 times to mix.

B. Streak Purifying Rehydrated, Lyophilized Bacteria

1. Label an agar plate with initials, date, and bacterial genus and species, and strain identifier (i.e.: *Arthrobacter globiformis*_NRRL_B-2979).
2. Transfer 10 μ L of the rehydrated material into the first quadrant of a streak plate.
3. Streak the 10 μ L for isolated colonies.
4. Place into incubator, checking for colonies roughly every 24 hours.
5. From the grown streak plate, pick a single colony and inoculate liquid growth media for a culture.
 - a. Monitor plate for contamination, colony morphology, and isolated colonies.
 - b. We recommend that you grow a minimum of 10 ml of culture to freeze.
 - c. These will be your frozen stock cultures.
6. Proceed to step E. to prepare your stock cultures.




Figure B: 1. Streak plate with single colonies. **2.** Freezer tube labeled with a freezer ready label, for a liquid culture freeze-down.

C. Check the culture for purity

1. On a non-antibiotic nutrient agar plate, streak out a loopful of your culture using aseptic technique.
2. Allow to grow for a week, checking each plate daily for purity (e.g.: single colony morphotype, growth rate).
 - a. For quick review, actinobacteria will not be fully grown in 24 hours. If you have colonies at 24 hours, suspect contamination and do not proceed.

NOTE:

- Because you do not want to delay making your frozen stock cultures, continue with the following steps at the same time. Do not wait for these results to continue.
- When initiating use of a stock culture, repeat a purity check before using any cultures that you may make. Always start a working culture from an isolated colony on a streak plate.

D. Preparing a Frozen Stocks of the Rehydrated Bacteria

Note: The generation number of the original material can be denoted as P.

1. Transfer the remaining rehydrated material and 40% glycerol (1:1) to a -80°C-compatible tube (hereafter freezer tube)
2. Prepare a label for the freezer tube (label) with the details of the bacteria and strain, as well as "Gen. P" and the date (Fig. B-2). Place this tube in the -80°C freezer, to be saved when all other stock cultures fail.
3. When preparing a streak plate from the frozen stock, do not allow the tube to thaw completely. Using a sterile loop, scrape the top of the freeze to obtain enough to streak out on a plate. Prepare the streak plate as normal.
4. We recommended preparing streak plates from -80°C in duplicate. Replace the tube in the -80°C as quickly as possible.

E. Preparing Multiple Frozen Stocks for Future Use

Note: The generation number of first generation cells can be denoted as F1.

1. Label a freezer tube (label) with the details of the bacteria and strain (Section B). Denote as F1 generation and add the date (Fig.B-2)
2. Add 500 µL of 40% glycerol to each tube.
3. Add 500 µL of your culture from step B-5 to each tube.
4. Seal the tube and place in -80°C.

Procedures: Phage

A. Lyophilized Phage Rehydration

1. With aseptic technique, transfer phage buffer into the vial.
2. The lyophilized cake will disintegrate almost immediately.
3. Gently pipette up and down 3-5 times to mix.

B. Plaque Purifying Rehydrated, Lyophilized Phage

1. Amplify your phage using the rehydrated sample. Follow protocol [Chapter 7: Phage amplification.](#)

C. Preparing a Frozen Stock of the Rehydrated Phage

1. With the remainder of the rehydrated phage - follow protocol [Chapter 7: Phage amplification.](#)

Created 5/7/2025
By E. Amarh/E. Ginser